

Immunophenotyping of tumor-specific CD8 T cells using high-dimensional mass cytometry

Michael Fehlings¹, Brian Abel², Alessandra Nardin¹, Evan Newell³, Mahesh Yadav⁴

¹immunoSCAPE, Singapore, ²immunoSCAPE, USA, ³Fred Hutchinson Cancer Research Center, USA, ⁴Genentech, Inc, USA

immunoSCAPE

Background: There is strong evidence that immunotherapy-mediated tumor rejection can be driven by tumor-specific CD8+ T cells reinvigorated to recognize neoantigens derived from tumor somatic mutations. Thus, the frequencies or characteristics of tumor-reactive, mutation-specific CD8+ T cells could be used as biomarkers of an anti-tumor response. However, such neoantigen-specific T cells are difficult to reliably identify due to their low frequency in peripheral blood and wide range of potential epitope specificities.


immunoSCAPE leverages the high-dimensional immune profiling capabilities of mass cytometry combined with a unique technology for the identification of antigen-specific T cells to address this challenge. By applying CyTOF in conjunction with combinatorial peptide-MHC tetramer staining and high-performance dimensional analysis tools, we can screen and identify rare neoantigen-specific T cells while concurrently performing deep phenotypic characterization.

Methods: Peripheral blood mononuclear cells (PBMC) from 14 non-small cell lung cancer (NSCLC) patients were collected pre- and post-treatment with the anti-PD-L1 antibody, atezolizumab, for characterization. Using whole exome and RNA sequencing we identified tumor neoantigens that are predicted to bind to major histocompatibility complex class-I (MHC-I) and utilized mass cytometry, together with cellular 'barcoding', to profile immune cells from patients with objective response to therapy (PR, n=8) and those with progressive disease (PD, n=6). In parallel, a highly-multiplexed combinatorial tetramer staining was used to screen antigen-specific CD8+ T cells in peripheral blood for 782 candidate tumor neoantigens and 71 known viral-derived control peptide epitopes across all patient samples.


Patient Response	PR	PD
Patients	8	6
Total NeoAgs predicted	487	302
Mean number of NeoAgs/patient	61	50

Summary of patient samples and characteristics (Fehlings *et al.*, JITC. 2019)


targetSCAPE: Large-scale epitope screening (100's) and deep antigen-specific T cell profiling *ex vivo*




No difference in bulk CD8+ T cells phenotype at baseline or following treatment between responders and non-responders



Neoantigen-specific T cells are enriched in patients responding to atezolizumab treatment




Neoantigen-specific T cells in responders have an effector phenotype similar to late-differentiated CMV-specific T cells



Conclusions:

- No significant treatment- or response-associated differences in bulk CD8+ T cell phenotype, but important to note that this is a small study
- Used an unbiased and sensitive method for *ex vivo* identification and characterization of rare antigen-specific T cells
- Most CD8+ T cell responses to neoAgs (13/20) are detected in patients with objective response to treatment, with 7 of 20 hits identified in patients with progressive disease
- In responders, most neoAg-specific T cells have an effector phenotype similar to late-differentiated CMV or effector-memory EBV-specific T cells
- The association between neoAg-specific T cell phenotype and treatment response should be confirmed in a larger cohort



Contact:
 Brian Abel, PhD, MBA
 immunoSCAPE, Inc, San Francisco
 Email: brian.abel@immunoscape.com
 Tel: (628) 999 3444

